

Product Datasheet

Phospho-CBL (Ser669) Antibody (orb1733105)

Description	This gene is a proto-oncogene that encodes a RING finger E3 ubiquitin ligase. The encoded protein is...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser669 of human CBL
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 100 kDa; Observed MW: 120 kDa
Uniprot ID	P22681
Entrez	867
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.