

Product Datasheet

Phospho-Smad2 (Ser255) Antibody (orb1733120)

Description	The protein encoded by this gene belongs to the SMAD, a family of proteins similar to the gene produ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser255 of human Smad2
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 52 kDa; Observed MW: 60 kDa
Uniprot ID	Q15796
Entrez	4087
Dilution Range	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.