

Product Datasheet

Phospho-GSK3 beta (Ser9) Antibody (orb1733123)

Description	Glycogen synthase kinase-3 (GSK3) is a proline-directed serine-threonine kinase that was initially i...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser9 of human GSK3 beta
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Uniprot ID	P49841
Entrez	2932
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.