

## Product Datasheet

### Phospho-RSK1 p90 (Thr359/Ser363) Antibody (orb1733124)

## Description

Rsk1 is a member of a family of

which includes Rsk1, Rsk2 and Rs...

## Species/Host

Rabbit

## Reactivity

Human, Mouse, Rat

## Conjugation

Unconjugated

## Tested

IP, WB

## Applications

## Immunogen

A synthetic phosphopeptide corresponding to residues surrounding Thr359/Ser363 of human RSK1 p90

## Preservatives

50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin

## Concentration

1mg/ml

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Note

For research use only

## Isotype

Rabbit IgG

## Clonality

Monoclonal

## MW

Calculated MW: 83 kDa; Observed MW: 90 kDa

## Uniprot ID

[Q15418](#)

## Entrez

[6195](#)

## Dilution Range

IP: 1/20; WB: 1/500-1/1000

## Expiration Date

12 months from date of receipt.