

## Product Datasheet

### Phospho-CBL (Tyr774) Antibody (orb1733125)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | This gene is a proto-oncogene that encodes a RING finger E3 ubiquitin ligase. The encoded protein is... |
| <b>Species/Host</b>        | Rabbit                                                                                                  |
| <b>Reactivity</b>          | Human, Rat                                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                                            |
| <b>Tested Applications</b> | IP, WB                                                                                                  |
| <b>Immunogen</b>           | A synthetic phosphopeptide corresponding to residues surrounding Tyr774 of human CBL                    |
| <b>Preservatives</b>       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin              |
| <b>Concentration</b>       | 1mg/ml                                                                                                  |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Isotype</b>             | Rabbit IgG                                                                                              |
| <b>Clonality</b>           | Monoclonal                                                                                              |
| <b>MW</b>                  | Calculated MW: 100 kDa; Observed MW: 120 kDa                                                            |
| <b>Uniprot ID</b>          | <b>P22681</b>                                                                                           |
| <b>Entrez</b>              | <b>867</b>                                                                                              |
| <b>Dilution Range</b>      | IP: 1/20; WB: 1/500-1/1000                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                         |