

Product Datasheet

Phospho-Presenilin 2 (Ser330) Antibody (orb1733130)

Description

Probable catalytic subunit of the

endoprotease complex that catalyzes th...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

A synthetic phosphopeptide corresponding to residues surrounding Ser330 of human Presenilin 2

Preservatives

50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin

Concentration

1mg/ml

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Monoclonal

MW

Calculated MW: 50 kDa; Observed MW: 23 kDa

Uniprot ID

P49810

Entrez

5664

Dilution Range

WB: 1/500-1/1000

Expiration Date

12 months from date of receipt.