

Product Datasheet

Phospho-Smad3 (Ser425) Antibody (orb1733143)

Description	Smad3 transcription factor phosphorylated and activated by TGF-beta-type receptors. A receptor-regul...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide conjugated to KLH.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P84022
Entrez	4088
Dilution Range	IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.