

## Product Datasheet

### COL1A1 Antibody (orb1733149)

|                            |                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | <span style="font-family:SimHei;">&lt;span style="font-family:"color:#000000;font-size:14px;"&gt;This gene...</span> |
| <b>Species/Host</b>        | Mouse                                                                                                                |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                                         |
| <b>Tested Applications</b> | IHC                                                                                                                  |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.03% Proclin 300, pH 7.3.                                  |
| <b>Concentration</b>       | 1mg/ml                                                                                                               |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                             |
| <b>Note</b>                | For research use only                                                                                                |
| <b>Isotype</b>             | Mouse IgG1                                                                                                           |
| <b>Clonality</b>           | Monoclonal                                                                                                           |
| <b>MW</b>                  | 139 kDa                                                                                                              |
| <b>Entrez</b>              | <b>1277</b>                                                                                                          |
| <b>Dilution Range</b>      | IHC: 1/100-1/200                                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                      |