

Product Datasheet

DiMethyl-Histone H3 (Lys27) Antibody (orb1733156)

Description	H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA acces...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Synthetic Peptide of Histone H3 (Di Methyl Lys27)
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Uniprot ID	P68431
Entrez	8350
Dilution Range	IP: 1/20; WB: 1/500-1/1000; ICC: 1/50- 1/200
Expiration Date	12 months from date of receipt.