

Product Datasheet

Cleaved-Caspase 3 p12 Antibody (orb1733170)

Description	Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Cas...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant protein of human Caspase-3 p12
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 32 kDa; Observed MW: 32,12 kDa
Uniprot ID	P42574
Entrez	836
Dilution Range	WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.