

## Product Datasheet

# Acetyl-Alpha Tubulin (Lys40) Antibody (orb1733176)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | TUBA1A Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchan... |
| <b>Species/Host</b>        | Mouse                                                                                                   |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                            |
| <b>Tested Applications</b> | IHC, WB                                                                                                 |
| <b>Immunogen</b>           | Synthetic peptide conjugated to KLH.                                                                    |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.                    |
| <b>Concentration</b>       | 1mg/ml                                                                                                  |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Isotype</b>             | Mouse IgG1                                                                                              |
| <b>Clonality</b>           | Monoclonal                                                                                              |
| <b>MW</b>                  | Calculated MW: 50 kDa; Observed MW: 50 kDa                                                              |
| <b>Uniprot ID</b>          | <b>P68363</b>                                                                                           |
| <b>Entrez</b>              | <b>10376</b>                                                                                            |
| <b>Dilution Range</b>      | WB: 1/500-1/1000; IHC: 1/50-1/100                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                         |