

Product Datasheet

EPO Antibody (orb1733199)

Description	Human erythropoietin is member of the EPO/TPO family and encodes a secreted, glycosylated cytokine h...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Purified recombinant fragment of human EPO expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	21kDa
Entrez	2056
Dilution Range	WB: 1/500 - 1/2000; ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.