


Product Datasheet

PRDM1/Blimp-1 Antibody (orb1733242)

Description	PRDM1/Blimp-1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 350 and 400
Target	PRDM1/Blimp-1
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	075626
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Detection of human PRDM1/Blimp-1 by west...