

Product Datasheet

Human ProS1 Protein, His Tag (orb1733275)

Description

This gene is a member of the leukocyte

which is found in ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

ProS1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 23.9 kDa. The protein migrates as 33-43 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human ProS1 Protein, His Tag (orb1733275) is expressed from human 293 cells (HEK293). It contains AA Leu 476 - Ser 676 (Accession # P07225).

MW

23.9 kDa

Uniprot ID

P07225

Expiration Date

6 months from date of receipt.