

Product Datasheet

Human TNF-beta / Lymphotoxin-alpha Protein, His Tag (MALS verified) (orb1733281)

Description

Human TNF-beta Protein, His Tag is

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TNF-beta

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 20.5 kDa. The protein migrates as 26-30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human TNF-beta Protein, His Tag (orb1733281) is expressed from human 293 cells (HEK293). It contains AA Leu 35 - Leu 205 (Accession # P01374).

MW


20.5 kDa

Uniprot ID

P01374

Expiration Date

6 months from date of receipt.



Immobilized
Human
TNFR1, Fc
Tag (Cat.
No. B01374)