

Product Datasheet

SARS-CoV-2 - nucleocapsid protein - 40.6 kDa (orb1733564)

Description	SARS-CoV-2 - coronavirus nucleocapsid protein
Form/Appearance	liquid (Supplied in Tris-HCl 50 mM, pH 8.0, Sodium Chlorid 60 mM, Glutathione 10 mM, Sarcosyl 0, 25%, 50% Glycerols)
Storage	store at -20 °C. avoid freeze/thaw cycles
Note	For research use only
Application notes	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not be excluded.
Clonality	Recombinant
MW	Theoretical MW: about 40.6 kDa
Expiration Date	6 months from date of receipt.