

Product Datasheet

SARS-CoV-2 - nucleocapsid protein - 48.3 kDa (orb1733565)

Description	SARS-CoV-2 - coronavirus nucleocapsid protein
Form/Appearance	liquid (Supplied in Sodium Phosphate 12.5 mM, pH 8.0, Tris-HCl 20 mM, Sodium Chlorid 300 mM, Imidasol 62.5 mM, Sarcosyl 0.5%)
Storage	store at 4 °C
Note	For research use only
Application notes	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not be excluded.
Clonality	Recombinant
MW	Theoretical MW: about 48.3 kDa
Expiration Date	6 months from date of receipt.