

## Product Datasheet

### KlenTaq (orb1733699)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | KlenTaq is a truncated version of Taq DNA Polymerase, lacking its first 280 amino acids. The enzyme ... |
| <b>Form/Appearance</b> | liquid                                                                                                  |
| <b>Concentration</b>   | 5 units/ $\mu$ l                                                                                        |
| <b>Storage</b>         | store at -20 °C. avoid freeze/thaw cycles                                                               |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Clonality</b>       | Recombinant                                                                                             |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                         |