

Product Datasheet

Fluorescent 20 bp DNA Ladder (orb1734248)

Description	Fluorescent 20 bp DNA Ladder consists of 10 DNA fragments ranging from 20 bp to 200 bp. The prestain...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, SYBR® Green and Orange G; Color: orange
Concentration	105 ng/μl
Storage	store at -20 °C. store dark. Short term storage (up to 1 month) at 4 °C possible
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.