

Product Datasheet

3-Azido-L-alanine HCl (orb1734915)

Description	(S)-2-Amino-3-azidopropanoic acid hydrochloride
Form/Appearance	powder; Color: white
Storage	store at 4 °C. store dry
Note	For research use only
Formula	$C_3H_6N_4O_2$
Purity	mass identification conforms (ESI-MS)
MW	Theoretical MW: 130.11 g/mol; Detected MW: 130.05 g/mol
CAS Number	105661-40-3
SMILES	<chem>OC(=O)[C@H](CN=[N+]=[N-])N</chem>
Expiration Date	12 months from date of receipt.