

Product Datasheet

KRT20 / CK20 / Cytokeratin 20 Antibody (clone KS 20.8) (orb1735016)

Description	KRT20 / CK20 / Cytokeratin 20 Antibody (clone KS 20.8)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Cytokeratins isolated from human micro dissected villi of human duodenal mucosa.
Preservatives	0.05% Sodium Azide
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Application notes	Further information: This antibody may be diluted to a titer of 1:25 to 1:50 in an ABC method. The final dilution should be optimized by the researcher. An incubation period of 30 minutes at RT is suggested. Formalin-fixed/paraffin-embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6, prior to immunostaining.
Isotype	IgG2a,k
Clonality	Monoclonal
Dilution Range	IHC-P (1:50)
Expiration Date	12 months from date of receipt.