

Product Datasheet

CCN5 Antibody (clone 1E11) (orb1735046)

Description	CCN5 Antibody (clone 1E11)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Recombinant protein
Preservatives	PBS, pH 7.4
Concentration	1 mg/ml
Storage	Store at -20°C or lower. Aliquot to avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Western Blot: Transfected lysate, Recombinant protein.
Isotype	IgG2a,k
Clonality	Monoclonal
Dilution Range	ELISA, IHC, IHC-P (5 µg/ml), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.