

Product Datasheet

APOB / Apolipoprotein B Antibody (orb1735047)

Description	APOB / Apolipoprotein B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ID, IHC, IHC-P, WB
Immunogen	Mouse very low density lipoprotein.
Preservatives	0.02% Sodium Azide
Storage	Short term: store at 4°C. Long term: store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Suitable for use in Radial Immunodiffusion (RID) and Western Blot. RID End-Point Titer: 1:15.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ID, IHC, IHC-P (1:100), WB
Expiration Date	12 months from date of receipt.