



Product Datasheet

S100B Antibody (orb1735083)


Description	S100B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide /S100 calcium binding protein B
Target	S100B
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500 - 1:1000 Immunofluorescence 1:100 - 1:400 Immunohistochemistry (Paraffin) 1:100 Immunohistochemistry (Frozen) 1:100 - 1:200
Expiration Date	12 months from date of receipt.



This blot was produced using a 15% SDS-P...



Immunohistochemical analysis of paraffin...



Immunohistochemical of frozen sections. ...