

Product Datasheet

Rat Kappa light chain Antibody (orb1735179)

Description	Rat Kappa light chain Antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	CS, ELISA, IA, IAC
Target	Rat Kappa light chain
Preservatives	PBS/0,1%NaN3. PBS/0,1%NaN3
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Uniprot ID	P01835
Expiration Date	12 months from date of receipt.