

Product Datasheet

Rat IgG2b heavy chain Antibody (HRP) (orb1735223)

Description	Rat IgG2b heavy chain Antibody (HRP)
Species/Host	Mouse
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA, ELISpot, IA, IAC
Target	Rat IgG2b heavy chain
Preservatives	PBS/50%glycerol. PBS/50%glycerol
Form/Appearance	HRP
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Uniprot ID	P20761
Expiration Date	12 months from date of receipt.