

Product Datasheet

Rat IgG2b heavy chain Antibody (Biotin) (orb1735227)

Description	Rat IgG2b heavy chain Antibody (Biotin)
Species/Host	Mouse
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA, ELISpot, IA, IAC
Target	Rat IgG2b heavy chain
Preservatives	PBS/0,1%NaN3/50%glycérol. PBS/0,1%NaN3/50%glycérol
Form/Appearance	Biotin
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Uniprot ID	P20761
Expiration Date	12 months from date of receipt.