

Product Datasheet

Rat IgG2b heavy chain Antibody (FITC) (orb1735234)

Description	Rat IgG2b heavy chain Antibody (FITC)
Species/Host	Mouse
Reactivity	Rat
Conjugation	FITC
Tested Applications	CS, ELISA, ELISpot, IA
Target	Rat IgG2b heavy chain
Preservatives	PBS/0,1%NaN3/50%glycérol. PBS/0,1%NaN3/50%glycérol
Form/Appearance	FITC
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Uniprot ID	P20761
Expiration Date	12 months from date of receipt.