

Product Datasheet

Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain Antibody (FITC) (orb1735295)

Description	Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain Antibody (FITC)
Species/Host	Mouse
Reactivity	Rat
Conjugation	FITC
Tested Applications	ELISA, IA
Target	Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain
Preservatives	PBS/0,1%NaN3/50%glycérol. PBS/0,1%NaN3/50%glycérol
Form/Appearance	FITC
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1,IgM,IgG2a kappa
Expiration Date	12 months from date of receipt.