

Product Datasheet

Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain Antibody (orb1735299)

Description	Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain Antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IA
Target	Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain
Preservatives	PBS/0,1%NaN3. PBS/0,1%NaN3
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1, IgG2a kappa
Expiration Date	12 months from date of receipt.