

## Product Datasheet

# Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain Antibody (orb1735301)

|                            |                                                  |
|----------------------------|--------------------------------------------------|
| <b>Description</b>         | Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain Antibody |
| <b>Species/Host</b>        | Mouse                                            |
| <b>Reactivity</b>          | Rat                                              |
| <b>Conjugation</b>         | Unconjugated                                     |
| <b>Tested Applications</b> | ELISA, IA                                        |
| <b>Target</b>              | Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain          |
| <b>Preservatives</b>       | PBS/0,1%NaN3.<br>PBS/0,1%NaN3                    |
| <b>Form/Appearance</b>     | Purified                                         |
| <b>Concentration</b>       | 1mg/mL                                           |
| <b>Storage</b>             | -20°C                                            |
| <b>Note</b>                | For research use only                            |
| <b>Isotype</b>             | IgG1, IgG2a kappa                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                  |