

## Product Datasheet

### Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain Antibody (Biotin) (orb1735306)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain Antibody (Biotin) |
| <b>Species/Host</b>        | Mouse                                                     |
| <b>Reactivity</b>          | Rat                                                       |
| <b>Conjugation</b>         | Biotin                                                    |
| <b>Tested Applications</b> | ELISA, IA                                                 |
| <b>Target</b>              | Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain                   |
| <b>Preservatives</b>       | PBS/0,1%NaN3/50%glycérol.<br>PBS/0,1%NaN3/50%glycérol     |
| <b>Form/Appearance</b>     | Biotin                                                    |
| <b>Concentration</b>       | 1mg/mL                                                    |
| <b>Storage</b>             | -20°C                                                     |
| <b>Note</b>                | For research use only                                     |
| <b>Isotype</b>             | IgG1, IgG2a kappa                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |