

Product Datasheet

Rat IgE heavy chain Antibody (orb1735310)

Description	Rat IgE heavy chain Antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	CS, ELISA, ELISpot, IA, IAC, WB
Target	Rat IgE heavy chain
Preservatives	PBS/0,1%NaN3. PBS/0,1%NaN3
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Uniprot ID	P01855
Expiration Date	12 months from date of receipt.