

## Product Datasheet

### Rat IgA heavy chain Antibody (orb1735340)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | Rat IgA heavy chain Antibody    |
| <b>Species/Host</b>        | Mouse                           |
| <b>Reactivity</b>          | Rat                             |
| <b>Conjugation</b>         | Unconjugated                    |
| <b>Tested Applications</b> | ELISA, ELISpot, IA, IAC, WB     |
| <b>Target</b>              | Rat IgA heavy chain             |
| <b>Preservatives</b>       | PBS/0,1%NaN3.<br>PBS/0,1%NaN3   |
| <b>Form/Appearance</b>     | Purified                        |
| <b>Concentration</b>       | 1mg/mL                          |
| <b>Storage</b>             | -20°C                           |
| <b>Note</b>                | For research use only           |
| <b>Isotype</b>             | IgG1 kappa                      |
| <b>Expiration Date</b>     | 12 months from date of receipt. |