

## Product Datasheet

### Rat kappa + lambda light chain Antibody (orb1735350)

|                            |                                         |
|----------------------------|-----------------------------------------|
| <b>Description</b>         | Rat kappa + lambda light chain Antibody |
| <b>Species/Host</b>        | Mouse                                   |
| <b>Reactivity</b>          | Rat                                     |
| <b>Conjugation</b>         | Unconjugated                            |
| <b>Tested Applications</b> | ELISA, IA                               |
| <b>Target</b>              | Rat kappa + lambda light chain          |
| <b>Preservatives</b>       | PBS/0,1%NaN3.<br>PBS/0,1%NaN3           |
| <b>Form/Appearance</b>     | Purified                                |
| <b>Concentration</b>       | 1mg/mL                                  |
| <b>Storage</b>             | -20°C                                   |
| <b>Note</b>                | For research use only                   |
| <b>Isotype</b>             | IgG1 kappa                              |
| <b>Uniprot ID</b>          | <b>P01681</b>                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.         |