

Product Datasheet

Rat kappa + lambda light chain Antibody (orb1735351)

Description	Rat kappa + lambda light chain Antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IA
Target	Rat kappa + lambda light chain
Preservatives	PBS/0,1%NaN ₃ . PBS/0,1%NaN ₃
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Uniprot ID	P01681
Expiration Date	12 months from date of receipt.