

Product Datasheet

Rat kappa + lambda light chain Antibody (HRP) (orb1735352)

Description	Rat kappa + lambda light chain Antibody (HRP)
Species/Host	Mouse
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA, IA
Target	Rat kappa + lambda light chain
Preservatives	PBS/50%glycerol. PBS/50%glycerol
Form/Appearance	HRP
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Uniprot ID	P01681
Expiration Date	12 months from date of receipt.