

Product Datasheet

Dinitrophenyl hapten Antibody (HRP) (orb1735372)

Description	Dinitrophenyl hapten Antibody (HRP)
Species/Host	Mouse
Conjugation	HRP
Tested Applications	ELISA, WB
Target	Dinitrophenyl hapten
Preservatives	PBS/50%glycerol. PBS/50%glycerol
Form/Appearance	HRP
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG3 kappa
Expiration Date	12 months from date of receipt.