

Product Datasheet

Dinitrophenyl hapten Antibody (Biotin) (orb1735396)

Description	Dinitrophenyl hapten Antibody (Biotin)
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Target	Dinitrophenyl hapten
Preservatives	PBS/0,1%NaN3/50%glycérol. PBS/0,1%NaN3/50%glycérol
Form/Appearance	Biotin
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG2a kappa
Expiration Date	12 months from date of receipt.