

Product Datasheet

Rat Kappa light chain of IgK-1b allotype Antibody (FITC) (orb1735414)

Description	Rat Kappa light chain of IgK-1b allotype Antibody (FITC)
Species/Host	Rat
Reactivity	Rat
Conjugation	FITC
Tested Applications	CS, IA, IAC
Target	Rat Kappa light chain of IgK-1b allotype
Preservatives	PBS/0,1%NaN3/50%glycérol. PBS/0,1%NaN3/50%glycérol
Form/Appearance	FITC
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Expiration Date	12 months from date of receipt.