

Product Datasheet

Rabbit IgG heavy chain Antibody (orb1735430)

Description	Rabbit IgG heavy chain Antibody
Species/Host	Rat
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	CS, ELISA, IA, IAC
Target	Rabbit IgG heavy chain
Preservatives	PBS/0,1%NaN ₃ . PBS/0,1%NaN ₃
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG2a kappa
Uniprot ID	P01870
Expiration Date	12 months from date of receipt.