

Product Datasheet

Mouse Kappa lighth chain Antibody (orb1735475)

Description	Mouse Kappa lighth chain Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IA, IAC
Target	Mouse Kappa lighth chain
Preservatives	PBS/0,1%NaN3. PBS/0,1%NaN3
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Uniprot ID	P01837
Expiration Date	12 months from date of receipt.