

Product Datasheet

Mouse Kappa ligth chain Antibody (HRP) (orb1735486)

Description	Mouse Kappa ligth chain Antibody (HRP)
Species/Host	Rat
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA, IA, IAC, IH
Target	Mouse Kappa ligth chain
Preservatives	PBS/50%glycerol. PBS/50%glycerol
Form/Appearance	HRP
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG2a kappa
Uniprot ID	P01837
Expiration Date	12 months from date of receipt.