

Product Datasheet

Mouse Kappa ligth chain Antibody (Biotin) (orb1735490)

Description	Mouse Kappa ligth chain Antibody (Biotin)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, IA, IAC, IH
Target	Mouse Kappa ligth chain
Preservatives	PBS/0,1%NaN3/50%glycérol. PBS/0,1%NaN3/50%glycérol
Form/Appearance	Biotin
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG2a kappa
Uniprot ID	P01837
Expiration Date	12 months from date of receipt.