

Product Datasheet

Mouse IgG1,IgG2a, IgG2b, IgG3 heavy chain Antibody (orb1735493)

Description	Mouse IgG1,IgG2a, IgG2b, IgG3 heavy chain Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IA
Target	Mouse IgG1,IgG2a, IgG2b, IgG3 heavy chain
Preservatives	PBS/0,1%NaN3. PBS/0,1%NaN3
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Expiration Date	12 months from date of receipt.