

Product Datasheet

Mouse IgG1 heavy chain Antibody (FITC) (orb1735589)

Description	Mouse IgG1 heavy chain Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	ELISA, IA, IAC
Target	Mouse IgG1 heavy chain
Preservatives	PBS/0,1%NaN3/50%glycérol. PBS/0,1%NaN3/50%glycérol
Form/Appearance	FITC
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Uniprot ID	P01868
Expiration Date	12 months from date of receipt.