

Product Datasheet

Mouse IgG1 heavy chain Antibody (HRP) (orb1735607)

Description	Mouse IgG1 heavy chain Antibody (HRP)
Species/Host	Rat
Reactivity	Mouse
Conjugation	HRP
Tested Applications	CS, ELISA, ELISpot, IA
Target	Mouse IgG1 heavy chain
Preservatives	PBS/50%glycerol. PBS/50%glycerol
Form/Appearance	HRP
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG1 kappa
Uniprot ID	P01868
Expiration Date	12 months from date of receipt.