

Product Datasheet

Mouse IgD heavy chain Antibody (orb1735645)

Description	Mouse IgD heavy chain Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ELISpot, IA, IAC
Target	Mouse IgD heavy chain
Preservatives	PBS/0,1%NaN3. PBS/0,1%NaN3
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG2a kappa
Uniprot ID	P01881
Expiration Date	12 months from date of receipt.