

## Product Datasheet

### Mouse IgA heavy chain Antibody (FITC) (orb1735658)

|                            |                                                       |
|----------------------------|-------------------------------------------------------|
| <b>Description</b>         | Mouse IgA heavy chain Antibody (FITC)                 |
| <b>Species/Host</b>        | Rat                                                   |
| <b>Reactivity</b>          | Mouse                                                 |
| <b>Conjugation</b>         | FITC                                                  |
| <b>Tested Applications</b> | ELISA, IA, IAC                                        |
| <b>Target</b>              | Mouse IgA heavy chain                                 |
| <b>Preservatives</b>       | PBS/0,1%NaN3/50%glycérol.<br>PBS/0,1%NaN3/50%glycérol |
| <b>Form/Appearance</b>     | FITC                                                  |
| <b>Concentration</b>       | 1mg/mL                                                |
| <b>Storage</b>             | -20°C                                                 |
| <b>Note</b>                | For research use only                                 |
| <b>Isotype</b>             | IgM kappa                                             |
| <b>Uniprot ID</b>          | <b>P01878</b>                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                       |