

Product Datasheet

Human IgG heavy chain Antibody (FITC) (orb1735801)

Description	Human IgG heavy chain Antibody (FITC)
Species/Host	Rat
Reactivity	Human
Conjugation	FITC
Tested Applications	ELISA, IA, IAC
Target	Human IgG heavy chain
Preservatives	PBS/0,1%NaN3/50%glycérol. PBS/0,1%NaN3/50%glycérol
Form/Appearance	FITC
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG2c kappa
Expiration Date	12 months from date of receipt.